

# Module 3

## Introduction to Wearables- Time Series

Cynthia Sandor

Edmond J Safra Assistant Professor  
UKRI Future Leader Fellow  
UK Dementia Research Institute Group Leader



# Announcements

- Next workshop: 11th November, CXGH 1.14 (Glenister Lecture Theatre), Charing Cross Campus (W6 8RP.)
- Today: Introduction to time-series data + Tutorial 3.
- **If you haven't finished Tutorial 2, please do so.**
- Tomorrow: Hackathon expectations overview.
- Tuesday: All tutorial solutions released.
- Reminder: Notebook → Google → TA or MT

Lecture theatres



CXGH 1.14 - Glenister Lecture Theatre

# What We'll Do Today

---

- **Lecture:** Short introduction to Wearables – Times Series Data


**Tutorial / Challenge —** Intro to **time-series data**: clean signals and extract features from **ECG** and **accelerometer** data.

- Work with real biomedical signals (ECG + accelerometer).
- Preprocessing: handle missing data, smoothing, filtering, normalization.
- ECG: R-peak detection, RR intervals, heart rate variability.
- Accelerometer: windowing, basic movement features.
- Feature extraction in time and frequency domains.



Cecilia  
Rodriguez

# Why should we focus on digital markers?



# Monitoring: What do we currently use?

Clinical visits: UPDRS

3.1 SPEECH	
<p><b>Instructions to examiner:</b> Listen to the patient's free-flowing speech and engage in conversation if necessary. Suggested topics: ask about the patient's work, hobbies, exercise, or how he got to the doctor's office. Evaluate volume, modulation (prosody) and clarity, including slurring, palilalia (repetition of syllables), and tachyphemia (rapid speech, running syllables together).</p> <p>0: Normal: No speech problems.</p> <p>1: Slight: Loss of modulation, diction, or volume, but still all words easy to understand.</p> <p>2: Mild: Loss of modulation, diction, or volume, with a few words unclear, but the overall sentences easy to follow.</p> <p>3: Moderate: Speech is difficult to understand to the point that some, but not most, sentences are poorly understood.</p> <p>4: Severe: Most speech is difficult to understand or unintelligible.</p>	
SCORE	<input type="text"/>
<hr/>	
3.2 FACIAL EXPRESSION	
<p><b>Instructions to examiner:</b> Observe the patient sitting at rest for 10 seconds, while talking. Observe eye-blink frequency, masked facies or loss of facial expression, and parting of lips.</p> <p>0: Normal: Normal facial expression.</p> <p>1: Slight: Minimal masked facies manifested only by decreased frequency of eye-blinking.</p> <p>2: Mild: In addition to decreased eye-blink frequency, masked facies as well, namely fewer movements around the mouth, spontaneous smiling, but lips not parted.</p> <p>3: Moderate: Masked facies with lips parted some of the time when the patient is at rest.</p> <p>4: Severe: Masked facies with lips parted most of the time when the patient is at rest.</p>	
	

Active virtual assessment

## Rune Labs scores FDA 510(k) to monitor Parkinson's symptoms on Apple Watch

The StrivePD app uses Apple's Movement Disorder API to track tremors and uncontrolled body movement.

By Emily Olsen | June 14, 2022 | 11:42 am

SHARE  20.6K



## NeuroRPM obtains FDA 510(k) for Parkinson's monitoring on Apple Watch


The AI-powered app allows for remote monitoring of common Parkinson's disease symptoms such as bradykinesia, tremor and dyskinesia via an Apple Watch.

By Jessica Hagen | March 24, 2023 | 11:18 am


SHARE  268



# StrivePD



# PKG



# Smartwatch data can predict, up to 7 years in advance, who is going to develop Parkinson's disease.



nature medicine

Explore content ▾

About the journal ▾


Publish with us ▾



[nature](#) > [nature medicine](#) > [articles](#) > [article](#)

Article | [Published: 03 July 2023](#)


## Wearable movement-tracking data identify Parkinson's disease years before clinical diagnosis





# Wearables Time-Series Data

6-dimensional inertial measurement unit, or IMU.






3-axis accelerometer



3-axis gyroscope





# From Raw Signals to Clinical Meaning: Tremor and Bradykinesia




# Revealing Hidden Rhythms in Wearable Sensor Data: Frequency Domain

(b)




# Revealing Hidden Rhythms in Wearable Sensor Data: Sampling Rate



$$\Delta t = 0.01sec = 10ms \rightarrow f_s = \frac{1}{\Delta t} = 100Hz \quad | \quad Hz = \frac{1}{sec}$$







# Aliasing — when a signal is sampled below its true frequency










# Power Spectral Density (PSD)

how the power of a signal is distributed across different frequencies


Time Domain vs Frequency Domain ( $T = 5$  s,  $fs = 100$  Hz)



# When Frequencies Change





# Spectrogram: Time-Frequency representation







# Filtering: How We Clean, Focus, and Reveal the True Signal

“Mathematical operations that attenuate the amplitude of a subset of frequencies in the signal.”  
→ “Attenuate” means **reduce strength or remove power** at those frequencies.



# High Pass Filtering

`highpass(Z, 0.5, 100, Steepness=0.95)`



**High-pass filters remove slow drifts and preserve fast dynamics.** This is essential in real-world data (e.g., tremor signals, accelerometer readings, EEG) to eliminate slow baseline fluctuations due to posture, movement, or sensor drift. Filtering helps us **focus on the meaningful activity frequencies** while cleaning up the noise or trend.

# Low Pass Filtering: A moving average


A **moving average** is a way to **smooth a signal** — to make it less noisy and easier to see the main trend.

A moving average replaces each sample by the **average of its neighbors**:



$$y[n] = \frac{1}{k} \sum_{i=0}^{k-1} x[n-i]$$

That operation is exactly the same as a **convolution** between:

- your signal  $x[n]$ , and
- a small kernel  $h[n] = \frac{1}{k}[1, 1, 1, \dots, 1]$ .



# Low Pass Filtering: A moving average






# Beyond IMU - Smartwatch




Multi-sensor device:  
Accelerometer  
Ambient temperature  
Gyroscope  
Humidity  
**Photoplethysmography (PPG)**

**verily**



# Beyond Smartwatch: Time-Serie Data very common in Neuroscience



# Let's Code



Group 1	Group 2	Group 3	Group 4	Group 5	Group 6	Group 7	Group 8
Philip Press	Viola Certani	Duru Okay	Yossra Semoukh	Harsimran Kaur	Jenna Lin	Mulreann Hogan	Licia Jing
Margarida Neves	Anya Kaur Haddon	Chloe Kam	Rebekah Boume	Luca Pastorello	Mustafa Gunaydin	Edona Bajrami	Zishan Lin
Olivia Pownall	Analys Chia	Wenbo Liao	Hattle Oliver	Luoyuan Zhang	Laura Raklec	Francesca Murley-Holme	Stephanie Sun
Charlotte Yu	Shobana Chandrashekhar	Yunyi Gao	Thilshany Kuganeswaran	Ellie Carre	Hanna Al tessald	Ruben Thilagarajah-Fernar	Neera Gahir
Chun Hei Leung	Katie Hay	Anas Saleem	Hanyue Pang	Veronika Shevchenko	Allsija Dabasinskaitė	Felix Varenne	Chuyi Zhang
Andrea Fan	Nina Jeffrey	Chi U Chau	Zehab Ben Hallim	Amina Bououdine	Ema Ferrá	Ruofan CAO	Lili Yassin
Adellina Shahata	Asma Abdullahi	YINUO Wang	Isabella Coloru	Krystal Tan	Temi Laina	Xinrui Fan	Keya Tanwani
Tanaka Udugama Jay Veer Bley		Mari Hronska	Lucas Yebra Garcia	Sarah Kurbanov			

[https://github.com/Sandoretal/Module\\_3/tree/main/tutorial\\_3](https://github.com/Sandoretal/Module_3/tree/main/tutorial_3)